



Original Article

Tissue and Plasma-Based Highly Sensitive Blocker Displacement Amplicon Nanopore Sequencing for *EGFR* Mutations in Lung Cancer

Patinya Akkhasutthikun^{1,2}, Pornchai Kaewsapsak^{1,2}, Pattaraporn Nimsamer^{1,3,4}, Pavit Klomkiew², Suthida Visedthorn², Pragwalai Chanchaem², Chinachote Teerapakpinyo⁵, Sunchai Payungporn^{1,2}, Sutima Luangdilok¹

¹Department of Biochemistry, Faculty of Medicine, Chulalongkorn University, Bangkok, ²Center of Excellence in Systems Microbiology, Faculty of Medicine, Chulalongkorn University, Bangkok, ³Division of Medical Bioinformatics, Research Department, Faculty of Medicine Siriraj Hospital, Mahidol University, Bangkok, ⁴Siriraj Long-Read Lab (Si-LoL), Faculty of Medicine Siriraj Hospital, Mahidol University, Bangkok, ⁵Chula GenePRO Center, Faculty of Medicine, Chulalongkorn University, Bangkok, Thailand

Purpose The epidermal growth factor receptor (*EGFR*) mutation is a widely prevalent oncogene driver in non-small cell lung cancer (NSCLC) in East Asia. The detection of *EGFR* mutations is a standard biomarker test performed routinely in patients with NSCLC for the selection of targeted therapy. Here, our objective was to develop a portable new technique for detecting *EGFR* (19Del, T790M, and L858R) mutations based on Nanopore sequencing.

Materials and Methods The assay employed a blocker displacement amplification (BDA)-based polymerase chain reaction (PCR) technique combined with Nanopore sequencing to detect *EGFR* mutations. Mutant and wild-type *EGFR* clones were generated from DNA from H1650 (19Del heterozygous) and H1975 (T790M and L858R heterozygous) lung cancer cell lines. Then, they were mixed to assess the performance of this technique for detecting low variant allele frequencies (VAFs). Subsequently, formalin-fixed, paraffin-embedded (FFPE) tissue and cell-free DNA (cfDNA) from patients with NSCLC were used for clinical validation.

Results The assay can detect low VAF at 0.5% mutant mixed in wild-type *EGFR*. Using FFPE DNA, the concordance rates of *EGFR* 19Del, T790M, and L858R mutations between our method and Cobas real-time PCR were 98.46%, 100%, and 100%, respectively. For cfDNA, the concordance rates of *EGFR* 19Del, T790M, and L858R mutations between our method and droplet digital PCR were 94.74%, 100%, and 100%, respectively.

Conclusion The BDA amplicon Nanopore sequencing is a highly accurate and sensitive method for the detection of *EGFR* mutations in clinical specimens.

Key words Non-small cell lung carcinoma, *EGFR* mutation, Blocker displacement amplification, Nanopore sequencing

Introduction

The epidermal growth factor receptor (*EGFR*) mutation is a common oncogene driver accounting for 50% of patients with non-small cell lung cancer (NSCLC) in East Asia [1]. Nowadays, *EGFR* tyrosine kinase inhibitors (*EGFR*-TKIs), namely erlotinib, gefitinib, and osimertinib, are used as the standard targeted therapy in NSCLC patients with *EGFR* mutations [2]. The 19Del and L858R *EGFR* mutations respond effectively to erlotinib and gefitinib treatment, inducing favorable outcomes [3,4]. However, the T790M mutation confers resistance to gefitinib or erlotinib therapy but remains sensitive to osimertinib, the third-generation *EGFR*-TKI [3,5]. Therefore, the detection of *EGFR* mutations is a standard routine testing to guide personalized treatment in clinical practice [2-4].

Cobas real-time polymerase chain reaction (PCR) is a stan-

dard method for detecting *EGFR* mutations in formalin-fixed, paraffin-embedded tissues (FFPE) with a sensitivity of at least 5% mutation levels, which is insufficient to detect cell-free DNA (cfDNA) containing low-level mutated DNA at less than 1% [6]. Instead, droplet digital PCR (ddPCR) is a technique that partitions the PCR reaction mixture into thousands of water-in-oil droplets, allowing for absolute quantification of target DNA with higher sensitivity than real-time PCR [7]. Therefore, ddPCR is used as the standard method for detecting *EGFR* mutations in cfDNA from blood [6,7]. These methods use an allele-specific forward primer designed to have a nucleotide at the 3' end perfectly matched with the mutant target [8]. However, this approach still has cross-reactivity with wild-type alleles [9]. To improve allele discrimination, the blocker methods, including peptide nucleic acid (PNA), locked nucleic acid (LNA), and blocker displace-

Correspondence: Sutima Luangdilok
Department of Biochemistry, Faculty of Medicine, Chulalongkorn University,
Bangkok 10330, Thailand
Tel: 66-22564482 Fax: 66-22564482 E-mail: sutima.l@chula.ac.th

Co-correspondence: Sunchai Payungporn
Center of Excellence in Systems Microbiology, Department of Biochemistry,
Faculty of Medicine, Chulalongkorn University, Bangkok 10330, Thailand
Tel: 66-22564482 Fax: 66-22564482 E-mail: sp.medbiochemcu@gmail.com

Received October 7, 2023 Accepted November 17, 2023
Published Online November 20, 2023

*Patinya Akkhasutthikun and Pornchai Kaewsapsak contributed equally to this work.

ment amplification (BDA), are used to avoid wild-type allele amplification [8,10].

PNA is a synthetic peptide polymer analog with nucleic acid bases as side chains, which can form hydrogen bonds to complementary DNA targets [11]. However, it lacks 3'OH, preventing DNA polymerase extension [11]. LNA is bridged nucleic acid that locks ribose moiety in the 3'-endo conformation, which prevents the cleavage mismatch nucleotide and enhances mismatch discrimination [12]. BDA uses a DNA oligonucleotide probe with internal three-carbon spacer, which can bind to wild-type DNA and prevent the binding of PCR primer [13,14]. Hence, BDA can limit the amplification of the wild-type alleles and effectively enrich low variant allele frequency (VAF) [9,13].

To identify the amplified product, the second-generation sequencing such as Illumina is commonly used [14]. However, these instruments are expensive and only available in centralized laboratories [14]. Recently Nanopore sequencing, a third-generation long-read sequencing, has gained more attention [14]. Its sequencer is portable and more affordable, allowing its usage in low resource-setting areas [15,16]. To combine Nanopore sequencing with PCR enrichment approach, BDA rather than PNA and LNA was chosen because it uses oligonucleotide to block the amplification of wild-type allele [13]. In contrast, LNA binds tightly to mutant alleles, resulting in amplification of mutant alleles [12,13]. However,

this amplicon contains LNA, which is not suitable for Nanopore sequencing [13]. For PNA, it is a peptide-based nucleotide analog, which is more difficult and expensive to synthesize than LNA and BDA probes [8,10].

To develop sensitive, specific, and cost-effective detection for *EGFR* 19Del, T790M, and L858R mutations, we combined BDA and Nanopore amplicon sequencing. We tested the performance of this assay in FFPE and cfDNA samples from NSCLC patients in the clinical setting (Fig. 1A and B).

Materials and Methods

1. Samples

The H1650 has 19Del mutation (E746_A750) in exon 19. The H1975 has T790M and L858R mutations in exons 20 and 21, respectively. Therefore, H1650 was used as a positive control for detecting *EGFR* 19Del, while H1975 was used for detecting *EGFR* T790M and L858R mutations. These cell lines were purchased from the American Type Culture Collection (ATCC, Manassas, VA) and cultured in Dulbecco's modified Eagle medium low glucose (Hyclone, Logan, UT) supplemented with 10% fetal bovine serum and 1% antibiotic-antimycotic (Gibco, Gaithersburg, MD) at 37°C in a humidified incubator with 5% CO₂.

The genomic DNA from H1650 and H1975 were extract-




Fig. 1. (A) The principle of blocker displacement amplification (BDA)-based polymerase chain reaction amplicon and Nanopore sequencing for epidermal growth factor receptor (*EGFR*) mutations. Formalin-fixed, paraffin-embedded (FFPE) and cell-free DNA (cfDNA) samples were performed based on BDA using blocker probe (light green line) bind with the wild-type allele (dark green circle) to prevent amplification. Conversely, it could not bind to the mutant allele (red circle), leading to amplification. Nanopore sequencing was used for following mutant allele enrichment. (B) Bioinformatic pipeline for variant detection (created with Biorender.com).

Table 1. PCR primers and blocker probes used for amplification of *EGFR* exons 19, 20, and 21

| Oligonucleotides name | Sequences (5'→3') | T _m (°C) | Product size (bp) |
|---------------------------|------------------------------------------------|---------------------|-------------------|
| <i>EGFR</i> -exon19-BDA-F | TTTCTGTTGGTGCTGATATTGCAGAAAGTTAAAATTCCTGTCGCTA | 61.3 | 180 |
| <i>EGFR</i> -exon19-BDA-R | ACTTGCCCTGTCGCTCTATCTTCGAAAAGGTGGGCCTGAGGTTTC | 66.1 | |
| <i>EGFR</i> -exon19-B | GTCGCTATCAAGGAATTAAGAGAAGCAA/iSpC3//iSpC3/CA | 66.1 | |
| <i>EGFR</i> -exon20-BDA-F | TTTCTGTTGGTGCTGATATTGCCACCTCCACCGTGCARCTCA | 64.8 | 188 |
| <i>EGFR</i> -exon20-BDA-R | ACTTGCCCTGTCGCTCTATCTTCTCCCTTCCCTGATTACCTTTGC | 65.1 | |
| <i>EGFR</i> -exon20-B | GTGCAGCTCATCACGCAGCTC/iSpC3//iSpC3/AT | 58.6 | |
| <i>EGFR</i> -exon21-BDA-F | TTTCTGTTGGTGCTGATATTGCGCAGCATGTCAAGATCACAGATT | 64.4 | 181 |
| <i>EGFR</i> -exon21-BDA-R | ACTTGCCCTGTCGCTCTATCTTCCAGCCTGGTCCCTGGTGTC | 67.3 | |
| <i>EGFR</i> -exon21-BD | GATCACAGATTTGGGCTGGCCAAA/iSpC3//iSpC3/CT | 67.2 | |

Underline bases represent the adapter nucleotides for PCR barcoding. "F" represents forward primers. "R" represents reverse primers. "B" represents blocker probes. "/iSpC3//iSpC3/" represents a C3 spacer modified nucleotide at the 3' end. BDA, blocker displacement amplification; *EGFR*, epidermal growth factor receptor; PCR, polymerase chain reaction.

ed by GenUP gDNA Kit (Biotechrabbit, Berlin, Germany), according to the manufacturer's instructions. The quantity of DNA was subsequently measured using the Nanophotometer C40 (Implen, Munchen, Germany).

A cohort of 65 FFPE DNA and 50 cfDNA samples from NSCLC patients with known *EGFR* mutations by Cobas real-time PCR and ddPCR, respectively, were obtained from the Chula GenePRO center, Faculty of Medicine, Chulalongkorn University under the approval of the Institutional Review Board of Faculty of Medicine, Chulalongkorn University, Bangkok, Thailand, IRB No. 0299/65 (S1 Table for the patient's age, sex, cancer types, and other specimen details).

2. Preparation of positive controls for detection

The *EGFR* gene in exons 19, 20, and 21 covering 19Del (position chr7:155728-155864), T790M (position chr7:162339-162483), and L858R (position chr7:172777-172914), respectively, were PCR-amplified using genomic DNA from the cancer cell lines. All amplicons were purified by QIAquick PCR Purification kit (Qiagen, Hilden, Germany), according to the manufacturer's instructions. Then, all purified samples were ligated to pGEM-T easy vector (Promega, Madison, WI). Ligation reaction, a total volume of 10 μ L contained 2.5 μ L of 5 \times rapid ligation buffer, 1 μ L of T4 DNA ligase, 1 μ L of 50 ng pGEM T easy (TA) vector, and 5.5 μ L of DNA insert. Subsequently, all ligated products were incubated at 22°C for 60 minutes and transformed into *Escherichia coli* DH5 α . The recombinant plasmids containing wild-type and mutant *EGFR* were extracted by PureDirex Plasmid mini-prep kit (BIO-HELIX Co., Ltd., Taipei, Taiwan), following the instructions provided by the manufacturer. The *EGFR* gene sequence in the recombinant plasmids was verified by Sanger sequencing, as shown in S2 Fig.

3. BDA assay

The three primer pairs and blocker probes were used to amplify *EGFR* exon 19 (position chr7:155728-155864), exon 20 (position chr7:162339-162483), and exon 21 (position chr7:172777-172914) covering each mutant position (Table 1, S3 Fig.). Forward primers were placed close to the mutant nucleotides, with 6 to 14 nucleotides overlapping with blocker probe sequences. Reverse primers were designed to prevent dimer formation between the primers and blockers. Subsequently, the 5' end of primers were tailed with adapter nucleotide for tagged with barcode for Nanopore sequencing. The blocker probes were designed to complement with wild-type allele. The 3' end of the blocker was modified into a 3-carbon spacer to prevent polymerase extension.

In the BDA-based PCR reaction, a total volume of 10 μ L contained 1 μ L of DNA template, 1 μ L of 1 \times Standard Taq Reaction Buffer (New England Biolabs, Hitchin, UK), 0.2 μ L of 200 μ M dNTPs mixed, 0.1 μ L of 10 μ M of each primer (U2Bio, Seoul, Korea), 1 μ L of 10 μ M blocker (U2Bio), 0.125 μ L of 100 U Taq DNA polymerase, and 6.475 μ L of DEPC treated water. The following PCR settings were used for amplification: initial denature at 95°C for 10 minutes, 35 cycles of amplification (95°C for 10 seconds, 58°C for 2 minutes, 68°C for 30 seconds) and final extension at 68°C for 10 minutes. Then, the amplified products were tagged with different barcodes.

For the PCR barcoding step, a total volume of 60 μ L contained 3 μ L of amplified products, 1 \times Standard Taq Reaction Buffer (New England Biolabs), 200 μ M dNTPs mixed, 0.2 μ M of each barcode primer from PCR Barcoding Expansion 1-96 (EXP-PBC096) kit (Oxford Nanopore Technologies, Oxford, UK), 0.5 U Taq DNA polymerase, and DEPC treated water. The PCR conditions included initial denaturation at 95°C for 10 minutes, 10 cycles of amplification (95°C for 10 seconds, 55°C for 1 minute, 68°C for 1 second), and final extension at




Fig. 2. The blocker displacement amplification (BDA) amplicon Nanopore sequencing was evaluated for variant allele frequency (VAF) enrichment (red bar), compared to without the blocker probes (blue bar). (A) The VAF enrichment of each mutant epidermal growth factor receptor (*EGFR*) from lung cancer cell lines. The detection of low VAF from mutant *EGFR* was mixed in wild-type *EGFR* with different ratios input from 5% to 0.5%. The VAF after BDA enrichment is shown in exon 19 (B), exon 20 (C), and exon 21 (D).

68°C for 5 minutes. After that, the amplified samples were purified by QIAquick PCR Purification kit (Qiagen), according to the manufacturer's instructions. Finally, the concentrations of purified products were measured by Qubit HS DNA Assay Kit (Invitrogen, Carlsbad, CA) and Qubit 4 fluorometer (Thermo Fisher Scientific, East Grinstead, UK).

4. DNA library preparation and sequencing

All samples were pooled equally to final amount of 1 µg DNA in distilled water. The adapters were ligated to the pooled samples using SQK-LSK114 kit (ONT, London, UK). The DNA library was loaded into a flow cell ver. R10.4.1 (FLO-MIN114) and sequenced by MinION platform. Approximately 10,000 and 30,000 reads per sample were sequenced for FFPE and cfDNA, respectively.

5. Bioinformatics

The raw read data in the FAST5 files were converted to the FASTQ files by Guppy basecaller software ver. 6.0.7 ([https://](https://github.com/oxfordnanopore/guppy)

github.com/asadprohdhan/GPU-accelerated-guppy-basecalling). To mitigate the Nanopore sequencing error background, we used super accuracy Phred quality score (Q score) ≥ 15 settings for base calling that includes an error rate of 3.16% or less. Thus, the sequencing used minimum read depth at 1,000 to ensure base accuracy. Then, the filtered reads were checked for quality control of sequencing data by MinION-QC (https://github.com/roblanf/minion_qc) before demultiplexed and adapter-trimmed by Porechop ver. 0.2.4. After that, total demultiplexed reads were rechecked in quality and length of sequenced read by FastQC (<https://www.bioinformatics.babraham.ac.uk/projects/fastqc/>) and then aligned with the human reference genome (GRCh38) using the Minimapp2 aligner (<https://github.com/lh3/minimapp2>) to generate SAM files. Then, the files were sorted and converted to BAM files by Samtools (<https://www.htslib.org/>). The Integrated Genome Viewer (IGV) program ver. 2.14.0 (<https://software.broadinstitute.org/software/igv/>) was utilized to visualize mutations and calculate VAF.




Fig. 3. The calculated cutoff for detecting epidermal growth factor receptor (*EGFR*) mutations in exon 19 (A), exon 20 (C), and exon 21 (E) from formalin-fixed, paraffin-embedded DNA. For cell-free DNA, the calculated cutoff for detecting *EGFR* mutations in exon 19 (B), exon 20 (D), and exon 21 (F). The red dots represent mutant mixed in wild-type alleles with ranging ratios of 5% to 0.1%. The blue dots represent the wild-type patient samples. LOD, limits of detection.

6. Identification of *EGFR* mutations status

The Clinical and Laboratory Standards Institute (CLSI) EP17-A guidelines were used to determine limits of blank (LOB) and limits of detection (LOD) [17,18]. In this study, LOB stands for the lowest VAF from sequencing error. The LOD stands for the lowest enriched VAF from sequencing.

Wild-type patients from FFPE (n=5) and cfDNA (n=9) samples were down-sampled to calculate sequencing error based on the %VAF. The enriched VAF studies were conducted

by diluting each mutant *EGFR* ranging from 5%, 2.5%, 1%, 0.5%, and 0.1%, which mimic low VAF in clinical specimens. Samples with %VAF higher than the determined cutoff were classified as mutants. Conversely, samples with %VAF lower than the cutoff were classified as wild type.

$$\text{Sequencing error} = \text{mean}_{\text{wild-type sample}} + 1.645 (\text{SD}_{\text{wild-type sample}}) \quad (1)$$

$$\text{Enriched VAF} = \text{sequencing error} + 1.645 (\text{SD}_{\text{low concentration sample}}) \quad (2)$$

Table 2. The performance for detecting *EGFR* mutations from FFPE DNA samples compared to Cobas real-time PCR

| | BDA amplicon Nanopore sequencing | | |
|----------------------|-------------------------------------|-------------|-------------|
| | 19Del | T790M | L858R |
| Sensitivity (%) | 100 (16/16) | 100 (10/10) | 100 (18/18) |
| Specificity (%) | 97.96 (48/49) | 100 (55/55) | 100 (47/47) |
| Concordance rate (%) | 98.46 (64/65) | 100 (65/65) | 100 (65/65) |

The data in parentheses are the numbers used to calculate percentages. BDA, blocker displacement amplification; *EGFR*, epidermal growth factor receptor; FFPE, formalin-fixed, paraffin-embedded; PCR, polymerase chain reaction.

7. Assessment of the performance in clinical specimens

To assess the performance of this technique, the results from FFPE DNA and cfDNA samples were then compared with the real-time PCR and ddPCR, respectively. The analytical sensitivity was calculated by $TP/(TP+FN)$. The specificity was calculated by $TN/(TP+FP)$, and concordance rate was calculated by $(TP+TN)/(TP+FP+FN+TN)$. TP, TN, FN, and FP indicate true-positive, true-negative, false-negative, and false-positive.

Results

1. Developing blocker displacement amplicon Nanopore sequencing to detect *EGFR* mutations

(1) Evaluation of the assay for detecting cancer cell lines

To develop a new technique for detecting *EGFR* mutations, BDA-based PCR amplicon Nanopore sequencing was designed and tested in lung cancer cell lines (H1650 and H1975) harboring *EGFR* mutations. These cell lines were used as the template in the PCR reaction with and without blocker to compare VAF enrichment. The amplified amplicons of exon 19, 20, and 21 (S4 Fig.) were sequenced by Nanopore sequencing. The result showed that the technique could increase VAF ranging from 93.38%, 92.10%, and 87.98% of *EGFR* 19Del, T790M, and L858R mutations, respectively, which were more than the average of 16.79% of the samples enriched by without blocker (Fig. 2A, S5 Table). Furthermore, all nucleotide mutations detected by the assay could be viewed on IGV, with a few sample instances presented in S6 Fig.

(2) Evaluation of the assay for detecting low VAF

To evaluate the assay for detecting low VAF in clinical specimens, the mutant samples input with different ratios of 5%, 2.5%, 1%, and 0.5% in a total DNA sample of 106 copies

Table 3. The performance for detecting *EGFR* mutations from cfDNA samples compared to ddPCR

| | BDA amplicon Nanopore sequencing | | |
|----------------------|-------------------------------------|-------------|-------------|
| | 19Del | T790M | L858R |
| Sensitivity (%) | 100 (22/22) | 100 (6/6) | 100 (10/10) |
| Specificity (%) | 87.5 (14/16) | 100 (32/32) | 100 (28/28) |
| Concordance rate (%) | 94.74 (36/38) | 100 (38/38) | 100 (38/38) |

The data in parentheses are the numbers used to calculate percentages. BDA, blocker displacement amplification; cfDNA, cell-free DNA; ddPCR, droplet digital polymerase chain reaction; *EGFR*, epidermal growth factor receptor.

were sequenced. The results indicated that BDA assay led to higher VAF enrichment than without blocker, the average of 6.80, 4.69, and 2.04-fold for exons 19, 20, and 21, respectively (Fig. 2B-D). Moreover, the detection limit was established at the lowest ratio of 0.5% enriched to 3.57%, 1.45%, and 0.83% for the mutations in exons 19, 20, and 21, respectively. However, the assay for detecting 0% mutant *EGFR* was found to have background noise from PCR amplification errors and a base-calling accuracy of 96.84% based on a Q score ≥ 15 (S7 Table).

2. Validating in clinical tissue specimens

(1) Calculated cutoff for detection

To distinguish between the sequencing error and low VAF, we calculated cutoff for detection following equations 1 and 2. Each *EGFR* mutant was serially diluted from 5% to 0.1% VAF and sequenced. The result found that exon 19 showed a mean enriched VAF of 5.81% (range, 0.73% to 11.13%), for exon 20 showed a mean enriched VAF of 2.99% (range, 0.09% to 7.04%), and for exon 21 showed a mean enriched VAF of 1.79% (range, 0.06% to 4.45%). The calculated sequencing error for wild-type samples indicated values of 1.22%, 0.49%, and 0.29% of exon 19, 20, and 21, respectively. Thus, the calculated cutoff for this technique was 3.09%, 1.44%, and 0.74% for 19Del, T790M, and L858R, respectively (Fig. 3A, C, and E).

(2) Evaluation of VAF in *EGFR* mutations

BDA amplicon Nanopore sequencing was applied for detecting *EGFR* mutations in FFPE DNA samples. The results showed that the VAF for 19Del ranged from 13.47% to 93.33%, T790M ranged from 6.41% to 49.46%, and L858R ranged from 17.10% to 81.78%, as shown in S8A Fig. In addition, 35 of 65 samples from NSCLC patients were detected with *EGFR* mutations, including 19Del (12 patients), L858R (13 patients), and double mutations as 19Del, T790M (5 pati-

ents) and T790M, L858R (5 patients). These results were in concordance with detection by the standard method based on Cobas real-time PCR, as shown in Table 2.

3. Validating in clinical blood specimens

(1) Calculated cutoff for detection

The validation of the technique for detecting *EGFR* mutations in cfDNA samples revealed the assay failed to amplify the *EGFR* gene using the same approach as for FFPE (data not shown). Subsequently, we optimized the BDA method to amplify *EGFR* mutations from cfDNA (S9 Table for the method). To determine a cutoff for detecting low VAF, wild-type cfDNA samples were used to calculate the VAF from sequencing error. The results found that the VAF was 3.06%, 0.23%, and 0.21% for exon 19, 20, and 21, respectively. The calculated cutoff was 4.94% for 19Del, 1.18% for T790M, and 0.66% for L858R mutations (Fig. 3B, D, and F).

(2) Evaluation of VAF in *EGFR* mutations

Thirty-eight of 50 samples were successfully amplified and sequenced. The VAF assessment found that the assay could enrich 19Del in a range from 13.25% to 66%, T790M in a range from 3.71% to 41.39%, and L858R in a range from 4.70% to 70.74% (S8B Fig.). Moreover, 29 of 38 cfDNA samples were detected *EGFR* mutations, including 19Del (15 patients), L858R (3 patients) and double mutations as 19Del and T790M (4 patients), T790M and L858R (2 patients), 19Del and L858R (5 patients). The result contained 12 false-positive samples for 19Del and one sample for L858R. Subsequently, these false-positive samples were validated through Barcode Taq (BT) sequencing. The results confirmed 10 of 12 cfDNA samples as true-positive for 19Del and one sample as true-positive for the L858R mutation. These results were consistent with the outcomes obtained from the ddPCR-based detection method, as shown in Table 3.

Discussion

In this study, we developed the technique to enable low VAF from *EGFR* mutations detection by BDA amplicon Nanopore sequencing for cost-effective detection. The assay demonstrated that using BDA significantly increased VAF. Despite mutant *EGFR* is normally presented at a very low level in NSCLC patients, this technique achieves detection of at least 0.5% VAF, compared to the NGS identify at least 5% VAF [19].

Furthermore, the assay was validated to detect *EGFR* mutations in clinical samples. Using multiple PCR cycles and more blocker concentration may lead to error detection of 19Del in cfDNA [20,21]. These were observed from

calculated VAF from sequencing error of wild-type patients, resulting in the cutoff higher than detecting *EGFR* mutations in FFPE.

For the testing with FFPE DNA, one sample was identified as a false-positive for the 19Del mutation. Subsequently, this sample was confirmed through real-time PCR, revealing it to be wild-type *EGFR*. Consequently, the detection of 19Del had a high concordance rate of 98.46% and 100% for T790M and L858R mutations.

For testing with cfDNA samples, 12 of 50 samples failed to amplify and were subsequently excluded, which might have low levels of cfDNA in the blood or degraded DNA. The assay found that T790M and L858R were 100% consistent with the ddPCR. Moreover, two samples were false-positive for 19Del mutation, resulting in a concordance rate of 94.74%. However, there are a few possible reasons for the discordance rate. Firstly, some NSCLC patients had *EGFR* exon 19 deletion of 10 bases (2243_2252_10del), hence Cobas *EGFR* probe for 19Del (15 bases deletion) would fail to capture. Secondly, some patients actually have rare double L858R and 19Del mutations [21,22]. According to this finding, 19Del and L858R mutations can potentially modify the molecular conformation of the *EGFR*-TK ATP binding site, particularly in the active domain [21,23]. This alteration may lead to reduced binding of *EGFR*-TKIs to the *EGFR*-TK domain, ultimately affecting the response to *EGFR*-TKIs [24]. However, the potential impact of concurrent mutations in exon 19 and 21 on clinical outcomes remains unclear [22]. Lastly, the BDA could enrich low VAF, suggesting that the assay might be more sensitive than the standard technique [14].

In the raw reads for sequencing, after base-calling using a Q score ≥ 15 results in a loss of approximately 36.21% of the total reads. However, a coverage depth of 1,000 reads is sufficient to detect gene mutation. Thus, for detecting three regions within one sample, a minimum raw read of 5,000 and reads per sample for FFPE DNA, and 7,500 reads per sample for cfDNA are recommended.

This technique could identify additional coincidental mutations within the range of amplicons, compared with the standard method that can only detect known mutations using specifically designed primers and probes. Moreover, the approach is portable and more affordable compared to the second-generation sequencers, permitting its usage in the low resource-setting areas. However, the turn-around time for the Cobas real-time PCR or ddPCR method was roughly 2 hours, significantly shorter than the turn-around time of this technique (approximately 3 days).

Nevertheless, this method had some limitations. First, the method is developed for detecting only common *EGFR* mutations (19del, T790M, and L858R), which does not cover other rare *EGFR* mutations (S10 Table). Second, this study

used the SQK-LSK114 kit for long-read sequencing. However, ONT plans to release a sequencing kit for short fragments, which might be appropriate for sequencing the DNA extracted from FFPE and cfDNA.

In conclusion, the development of BDA Nanopore sequencing proves to be highly sensitive for detecting *EGFR* 19Del, T790M, and L858R mutations in NSCLC and could be applied in clinical settings.

Electronic Supplementary Material

Supplementary materials are available at Cancer Research and Treatment website (<https://www.e-crt.org>).

Ethical Statement

This study used FFPE DNA and cfDNA from leftover specimens from routine clinical care that were not individually identifiable. All patients provided informed consent for specimen collection and *EGFR* testing. The need for the informed consent of the participants for reviewing medical records was waived on the condition that the research involves no more than minimal risk to the patients and their privacy. The study was eligible for expedited review and approved by the International Review Board of the Faculty of Medicine, Chulalongkorn University (IRB No. 0299/65).

Author Contributions

Conceived and designed the analysis: Akkhasutthikun P, Kaewsapsak P, Payungporn S, Luangdilok S.

Collected the data: Akkhasutthikun P, Teerapakpinyo C, Luangdilok S.

Contributed data or analysis tools: Akkhasutthikun P.


Performed the analysis: Klomkiew P, Visedthorn S.


Wrote the paper: Akkhasutthikun P, Kaewsapsak P, Payungporn S, Luangdilok S.


Contributed data: Nimsamer P, Chanchaem P, Teerapakpinyo C, Luangdilok S.

ORCID iDs

Patinya Akkhasutthikun  : <https://orcid.org/0009-0001-7416-6156>

Pornchai Kaewsapsak  : <https://orcid.org/0000-0001-7921-0868>

Sunchai Payungporn  : <https://orcid.org/0000-0003-2668-110X>

Sutima Luangdilok  : <https://orcid.org/0000-0003-1415-2331>

Conflicts of Interest

Conflict of interest relevant to this article was not reported.

Acknowledgments

The research was funded by H.M. the King Bhumibhol Adulyadej's 72nd Birthday Anniversary Scholarship of Chulalongkorn University and Ratchadapiseksompotch Fund, Faculty of Medicine, Chulalongkorn University, Thailand [Grant number RA65/005 and GA66/032].

References

- Hirsch FR, Varella-Garcia M, Bunn PA Jr, Di Maria MV, Veve R, Bremmes RM, et al. Epidermal growth factor receptor in non-small-cell lung carcinomas: correlation between gene copy number and protein expression and impact on prognosis. *J Clin Oncol.* 2003;21:3798-807.
- Miller VA, Kris MG, Shah N, Patel J, Azzoli C, Gomez J, et al. Bronchioloalveolar pathologic subtype and smoking history predict sensitivity to gefitinib in advanced non-small-cell lung cancer. *J Clin Oncol.* 2004;22:1103-9.
- Gupta R, Dastane AM, Forozan F, Riley-Portuguez A, Chung F, Lopategui J, et al. Evaluation of *EGFR* abnormalities in patients with pulmonary adenocarcinoma: the need to test neoplasms with more than one method. *Mod Pathol.* 2009;22:128-33.
- Riely GJ, Pao W, Pham D, Li AR, Rizvi N, Venkatraman ES, et al. Clinical course of patients with non-small cell lung cancer and epidermal growth factor receptor exon 19 and exon 21 mutations treated with gefitinib or erlotinib. *Clin Cancer Res.* 2006;12:839-44.
- Bethune G, Bethune D, Ridgway N, Xu Z. Epidermal growth factor receptor (*EGFR*) in lung cancer: an overview and update. *J Thorac Dis.* 2010;2:48-51.
- Deepak S, Kottapalli K, Rakwal R, Oros G, Rangappa K, Iwahashi H, et al. Real-time PCR: revolutionizing detection and expression analysis of genes. *Curr Genomics.* 2007;8:234-51.
- Pinheiro LB, Coleman VA, Hindson CM, Herrmann J, Hindson BJ, Bhat S, et al. Evaluation of a droplet digital polymerase chain reaction format for DNA copy number quantification. *Anal Chem.* 2012;84:1003-11.
- Guha M, Castellanos-Rizaldos E, Makrigiorgos GM. DISSECT method using PNA-LNA clamp improves detection of *EGFR* T790m mutation. *PLoS One.* 2013;8:e67782.
- Yao J, Zhang Z, Huang X, Guo Y. Blocker displacement amplification mediated PCR based screen-printed carbon electrode biosensor and lateral flow strip strategy for CYP2C19*2 genotyping. *Biosens Bioelectron.* 2022;207:114138.
- Smith GD, Zhou L, Rowe LR, Jarboe EA, Collins BT, Bentz JS, et al. Allele-specific PCR with competitive probe blocking for sensitive and specific detection of *BRAF* V600E in thyroid fine-needle aspiration specimens. *Acta Cytol.* 2011;55:576-83.
- Pellestor F, Paulasova P. The peptide nucleic acids (PNAs), powerful tools for molecular genetics and cytogenetics. *Eur J*

- Hum Genet. 2004;12:694-700.
12. Latorra D, Campbell K, Wolter A, Hurley JM. Enhanced allele-specific PCR discrimination in SNP genotyping using 3' locked nucleic acid (LNA) primers. *Hum Mutat.* 2003;22:79-85.
 13. Xu J, Pu Y, Lin R, Xiao S, Fu Y, Wang T. PEAC: an ultrasensitive and cost-effective MRD detection system in non-small cell lung cancer using plasma specimen. *Front Med (Lausanne).* 2022;9:822200.
 14. Thirunavukarasu D, Cheng LY, Song P, Chen SX, Borad MJ, Kwong L, et al. Oncogene concatenated enriched amplicon Nanopore sequencing for rapid, accurate, and affordable somatic mutation detection. *Genome Biol.* 2021;22:227.
 15. Wang Y, Zhao Y, Bollas A, Wang Y, Au KF. Nanopore sequencing technology, bioinformatics and applications. *Nat Biotechnol.* 2021;39:1348-65.
 16. Norris AL, Workman RE, Fan Y, Eshleman JR, Timp W. Nanopore sequencing detects structural variants in cancer. *Cancer Biol Ther.* 2016;17:246-53.
 17. Mimoso ML, Al-Ameri W, Simpson JT, Nakhla M, Boissinot K, Munoz DG, et al. A novel approach to detect IDH point mutations in gliomas using Nanopore sequencing: test validation for the clinical laboratory. *J Mol Diagn.* 2023;25:133-42.
 18. Pierson-Perry JF, Vaks JE, Vore TE, Durham AP, Fischer C, Gutenbrunner C, et al. Evaluation of detection capability for clinical laboratory measurement procedures; approved guideline (EP17-A2). Wayne, PA: Clinical Laboratory Standards Institute; 2012.
 19. Lopez-Rios F, Angulo B, Gomez B, Mair D, Martinez R, Conde E, et al. Comparison of molecular testing methods for the detection of EGFR mutations in formalin-fixed paraffin-embedded tissue specimens of non-small cell lung cancer. *J Clin Pathol.* 2013;66:381-5.
 20. Keschull JM, Zador AM. Sources of PCR-induced distortions in high-throughput sequencing data sets. *Nucleic Acids Res.* 2015;43:e143.
 21. Malapelle U, Sirera R, Jantus-Lewintre E, Reclusa P, Calabuig-Farinas S, Blasco A, et al. Profile of the Roche cobas(R) EGFR mutation test v2 for non-small cell lung cancer. *Expert Rev Mol Diagn.* 2017;17:209-15.
 22. Wei Z, An T, Wang Z, Chen K, Bai H, Zhu G, et al. Patients harboring epidermal growth factor receptor (EGFR) double mutations had a lower objective response rate than those with a single mutation in non-small cell lung cancer when treated with EGFR-tyrosine kinase inhibitors. *Thorac Cancer.* 2014;5:126-32.
 23. Harrison PT, Vyse S, Huang PH. Rare epidermal growth factor receptor (EGFR) mutations in non-small cell lung cancer. *Semin Cancer Biol.* 2020;61:167-79.
 24. Liang H, Li C, Zhao Y, Zhao S, Huang J, Cai X, et al. Concomitant mutations in EGFR 19Del/L858R mutation and their association with response to EGFR-TKIs in NSCLC patients. *Cancer Manag Res.* 2020;12:8653-62.